


Fluorescence Microscopy

Instruments, Methods, Applications

by Karl-Friedrich Koch
Ernst Leitz GmbH Wetzlar



Preface

This booklet attempts to present the use of fluorescence microscopes in a concise and lucid form according to the latest developments of the instruments in this field, without however, going into a detailed explanation of theoretical and practical principles of fluorescence microscopy.

On the other hand, a comprehensive description is offered of the instruments available and their applications to the various methods of investigation based on fluorescence microscopy from visual observation to fluorescence photomicrography.

The booklet has therefore been written mainly for the user of LEITZ fluorescence microscopes who is already versed in the relevant techniques of preparation and would like to obtain specialized information about the use of fluorescence microscopic devices in his field of research.

A list of references supplements the technical description.

I wish to thank Mrs. Ursula Daun and Mr. Helmut Rühl for their assistance with my investigations of the excitation and fluorescence maxima of the fluorochromes and with the plotting of the filter curves; my thanks are also due to Dr. F. Walter, Head of the Laboratory for Applied Microscopy, for his unstinting help with the production of this brochure.

Karl-Friedrich Koch,
Applied Microscopy Laboratory,
Ernst Leitz GmbH.


Wetzlar, February 1972

Contents

| | | |
|-----|---|----|
| I | What is Fluorescence? | 4 |
| A | Primary fluorescence | 5 |
| B | Secondary fluorescence through fluorochroming | 5 |
| C | Immunofluorescence | 5 |
| D | Histochemical and cytochemical fluorescent staining | 5 |
| II | Microscopic instruments of fluorescence | |
| | microscopy | 6 |
| A | Light sources | 8 |
| | Technical data | 9 |
| | The use of light sources in colour photography | 10 |
| B | Filters | 10 |
| | Broad-band excitation | 10 |
| | Narrow-band excitation | 11 |
| | Selective excitation | 12 |
| | Filter curves | 13 |
| | Fluorescent dyes | 17 |
| | Filter combinations | 20 |
| C | Methods of illumination in fluorescence | |
| | microscopy | 21 |
| | Incident-light fluorescence | 21 |
| | Transmitted-light fluorescence | 23 |
| | Simultaneous methods | 23 |
| D | Objectives | 36 |
| E | Eyepieces | 37 |
| III | Fluorescence photomicrography | 38 |
| A | The specimen | 38 |
| B | The microscope | 38 |
| C | The camera | 38 |
| D | The photographic material | 39 |
| IV | References | 40 |
| V | List of LEITZ literature on the instruments mentioned in this booklet | 43 |
| VI | Index | 44 |

I What is Fluorescence?


Certain substances have the property, when irradiated with ultra-violet, violet, blue, or green light, of emitting radiation of their own whose wave length is longer than that of the exciting light (Stoke's Law). This emission phenomenon is called **luminescence**:



27627 a-512 R

We must however, distinguish between two forms of luminescence:

If after the exciting light is switched off the secondary emission persists for some time at appreciable intensity, we speak of **phosphorescence**:




Phosphorescence during
irradiation




Phosphorescence after the
light source has been
switched off

27626-512 R

If, however, the excited radiation persists only so long as the exciting light acts on the specimen, we speak of **fluorescence**:



Fluorescence during
irradiation only



No fluorescence after the
light source has been
switched off

27627-512 R

Key to symbols:

E = Exciting filter


Er = Exciting light

Em = Emitted light

L = Lamp housing with light source

A Primary fluorescence

A number of unstained specimens have the property of emitting fluorescent light when they are excited with short-wave light. This fluorescence is called **primary** or **auto-fluorescence**:



B Secondary fluorescence through fluorochroming

Many specimens, however, do not exhibit any primary fluorescence. If they are stained with certain dyes — the fluorochromes — tissues that do not appear, but are in fact, different, can be shown in contrasting colours and distinguished by fluorescence.

The light emission caused by fluorochromes is called **secondary fluorescence**.

The dyes acridine orange, auramine, coriphosphin, rhodamine, phosphin 3 R are some of the best known fluorochromes.

Compounds such as the antibiotic tetracycline, which is stored by certain tissues of the organism (bone, teeth) and has fluorescent properties, are fluorochromes in the wider sense.

C Immunofluorescence

The method of immunofluorescence is highly specific because of its genuine serological reactions. Compounds


such as fluorescein isothiocyanate (FITC) and tetramethyl rhodamine isothiocyanate (TRITC) are used as fluorescent labelling substances. These are coupled with the antibodies by co-valent bonds. They are made to react with the corresponding antigen (e.g. pathogenic agent). The antigen-antibody complex produced can be recognized by its fluorescence when it is irradiated with short-wave light.

D Histochemical and cytochemical fluorescent staining

Another method of preparation involves histochemical or cytochemical fluorescent staining. Fluorescent end products are obtained by treatment of the material to be examined with chemical reagents.

The best known method is the process of condensation of biogenic amines (catecholamine, hydroxytryptamine) with formaldehyde, developed by ERÄNKÖ (1955) and FALCK (1962).

The technique developed by CASPERSSON et al. (1969, 1970) of cytochemical staining of chromosomes with a mustard gas compound for the demonstration and measurement of DNA distribution also belongs to this category.



Untreated specimen
exhibits no fluorescence

Fluorescence after treatment
of the specimen (fluorochroming,
serological reaction, etc.)

II Microscopic instruments for fluorescence microscopy

Modern fluorescence microscopy relies on the use not only of specific methods of preparation and staining but also of a special fluorescence microscope. In addition the interest of research has recently been focused more and more on fluorometry as a quantitative method through measurement of the intensity of fluorescent radiation.

From these aspects the following stands have been developed for fluorescence microscopy:


The microscope for routine investigations in diagnostics, preferably for transmitted-light fluorescence microscopy.

The specialized fluorescence microscope for modern incident-light fluorescence excitation and investigation in transmitted light.

The microscope fluorometer, consisting of a research microscope and the MPV microscope photometer as basic instruments, as well as recording and evaluation devices.

The basic stands for fluorescence microscopy are the LEITZ DIALUX, ORTHOLUX, and ORTHOPLAN microscopes.

Microscope outfits for fluorometry are described in Section II/C of this booklet.



27646-512 R

Basic construction of the fluorescence microscope

- | | |
|-------------------|----------------------|
| 1 Light source | 4 Objective |
| 2 Exciting filter | 5 Suppression filter |
| 3 Condenser | 6 Eyepiece |

The exciting light is produced by a suitable light source. Appropriate exciting filters (primary filters) are inserted in the lamp housing; they transmit only the exciting radiation, which, through a condenser (transmitted-light excitation) or through the fluorescence vertical illuminator (incident-light excitation) reaches the specimen, where it excites fluorescence.

For the observation of the fluorescent radiation emitted by the specimen a suppression filter (secondary filter) matched for the wave length of the fluorescent light must be inserted in the beam. It absorbs (glass filter) or reflects (interference filter) the exciting radiation liable to damage the eye and produces a dark image background.

To concentrate the highest possible amount of exciting energy on the specimen, light source and condenser must be accurately centred. The instructions "Lamp Housing 250 (514-72)" "Lamp Housing 100/100 Z (514-119)" and "Darkground Condensers D 0.80/D 1.20 (513-31)" contain further information.


For optimum yield of fluorescence the combination of exciting and suppression filters matched for the specimen under investigation is also important.

Intensity comparison of various light sources suitable for fluorescence microscopy (referred to 6 v 30 W low-voltage lamp)


A Light sources

| Lamp | Use in incident- light fluorescence | Use in transmitted- light fluorescence | Fluoro- metry |
|--|--|---|------------------|
| 12 v 100 W tungsten halogen lamp | x | x (routine, FITC) | x |
| 75 W high pressure xenon lamp | x | -*) | x |
| 150 W high pressure xenon lamp | x | x | |
| 450 W high pressure xenon lamp | x | x | - |
| 50 W ultra-high pressure mercury lamp | x | x (routine) | - |
| 100 W ultra-high pressure mercury lamp | x | -*) | x |
| 200 W ultra-high pressure mercury lamp | x | x | - |

* these lamps are not suitable for transmitted-light excitation, since they do not fill the entrance pupil of the darkground condensers.



Spectral energy distribution of the 200 W ultra-high pressure mercury lamp as example of a suitable fluorescence light source



1 nm = 10 A. U. = 10^{-9} m.

Technical data of lamps suitable for fluorescence microscopy

Low-voltage lamps

| Lamp | Electrical power (W) | Voltage (V) | Lamp current (A) | Light flux (lm) | Colour temperature (°K) | Life in hours |
|----------------------------------|----------------------|-------------|------------------|-----------------|-------------------------|---------------|
| 12 v 100 W tungsten halogen lamp | 100 | 12 | 8 | 2800 | 3100 at 7.4 amps | ca 50 |

Xenon high-pressure lamps

| Lamp | Power (W) | Voltage (V) | Lamp current (A) | Light flux (lm) | Mean luminous density (sb) | Dimensions of the illuminated field (mmLxB) | Colour temperature (°K) | Life in hours |
|-------|-----------|-------------|------------------|-----------------|----------------------------|---|-------------------------|---------------|
| 75 W | 75 | 14 | 5,4 | 1 000 | 40 000 | 0,5 x 0,25 | 6000 | ca. 400 |
| 150 W | 150 | 20 | 7,5 | 3 000 | 15 000 | 2,2 x 0,5 | 6000 | ca. 1200 |
| 450 W | 450 | 18 | 25 | 13 000 | 35 000 | 2,7 x 0,9 | 6000 | ca. 2000 |

Ultra-high pressure mercury lamps

| Lamp | Power (W) | Voltage (V) | Lamp current (A) | Light flux (lm) | Mean luminous density (sb) | Dimensions of the illuminated field (mmLxB) | Life in hours |
|-------|-----------|------------------------|--------------------|-----------------|----------------------------|---|---------------|
| 50 W | 50 | L1 39..45 L2 34..39 | 1,3 1,45 | 2000 | 30 000 | 1,0 x 0,3 | ca. 100 |
| 100 W | 100 | 20 | 5 | 2200 | 170 000 | 0,25 x 0,25 | ca. 200 |
| 200 W | 200 | L1 61 ± 4 L2 53 ± 4 | 3,4–3,9 3,9–4,4 | 9500 | 33 000 | 2,2 x 0,5 | ca. 200 |

Light sources in colour photography

Fluorescence colour photography: All light sources are suitable (daylight colour film)

Other microscopic methods: for artificial-light colour film 12v 100W tungsten halogen lamp

For daylight colour film
high-pressure xenon lamps
75 W
150 W
450 W

Except in fluorescence colour photography the ultra-high pressure mercury lamps are not suitable for use with colour film.

B Filters

Not only the choice of a suitable light source, but also that of the filters and their combination depends on the nature of the specimen. Both must be matched for the fluorescence staining method used. This requires knowledge of the spectral values for the excitation and fluorescence (emission).


We generally distinguish between filter combinations for broad-band, narrow-band, and selective excitation.

Broad-band excitation

Broad-band fluorescence constitutes the classical method. Here the exciting light is filtered out with the aid of dyed glass filters (UG filters, BG filters) of great half-value width and the fluorescent light observed from a certain wavelength onwards through suppression filters.

The method is used in transmitted light as well as in incident light wherever the excitation and emission maxima of the fluorochrome are widely separated. With the aid of the broad transmission ranges of the exciting filters the fluorescence can be excited with sufficient energy in or near the maximum without the emitted fluorescent light being influenced by the exciting radiation. This means that the observer can use a suppression filter which, because of the position of the edge, practically completely suppresses the exciting radiation, but fully transmits the maximum of the fluorescent radiation and therefore the "true" colour of the fluorescent light.


The illustration shows the example of a typical case of acriflavine fluorescence.




The maxima of excitation and emission are separated by about 100nm. Thus a BG 12 (32mm) exciting filter and K 510 or K 530 suppression filter can be chosen as combination for effective broad-band fluorescence.

Narrow-band excitation

Narrow-band excitation is used wherever the excitation and emission maxima of fluorescence are closely spaced. Here the exciting light is filtered out with interference band filters of small half-value width if necessary combined with conventional glass filters. This permits the excitation of fluorescence in a narrow range around the excitation maximum, without disturbance of the fluorescence emission by exciting light.



The usual edge filters are used as suppression filters; for the strict selection of fluorescence within the range of the emission maximum interference selection filters are used. The illustration above shows an example of narrow-band fluorescence. It illustrates the situation during excitation of red fluorescence of tetra-methyl rhodamine isothiocyanate (TRITC) by green light of 546nm. The TRITC maxima are so closely spaced here that the primary filter combination of the S 546 interference band filter (type AL) with the BG 38 and BG 36 glass filters becomes necessary to cut off the exciting radiation within the range of the emission maximum of TRITC. The graph further demonstrates that, because of the position of the maxima and the properties of the primary filter combination it is not the yellow emission maximum, but an intense red-orange or red fluorescence (PLOEM 1969) that can be



monofluorescence is the most widely used at the moment. With the KP 490 interference filter the excitation is produced directly within the range of the excitation maximum of FITC (about 490nm). For the observation of the characteristic apple green FITC fluorescence the K 510 edge filter, and particularly the S 525 selection filter is used for suppression. The additional BG 38 colour glass filter prevents the low transmission of the KP filter in the red from becoming visible.


observed, since a K 590 or K 610 suppression filter must be inserted.

Selective excitation


Similar to narrow-band excitation, selective excitation is used wherever excitation and emission maxima of the fluorescence are closely spaced. It differs from narrow-band excitation in that asymmetrical interference filters with a particularly steep cut-off are used. These filters have their maximum transmission in the excitation maximum of the fluorescent dye. Immediately next to the maximum their transmission is strictly limited in the direction of the longer wave length by the cut-off. It is combined on the observer's side with the usual edge suppression filters or with interference selection filters. Selective excitation of FITC im-

Filter curves

Filters for UV excitation




1mm UGI UV exciting filter




2mm UG 1,
2+2mm UG 1


Filters for violet excitation



3mm BG 3 violet exciting filter




3+3mm BG 3 violet exciting filter




S 405 interference band filter


Filter for blue excitation




1.5mm BG 12 blue exciting filter



3mm BG 12 blue exciting filter




S 470 Interference band filter




KP 490 FITC filter

Filter for green excitation




2mm BG 36+S 546 green exciting filter

Red suppression filter




4mm BG 38 red suppression filter

Heat filter




2mm KG 1 heat absorption filter


Suppression filter




K 430, K 460



Suppression filter



Filter for fluorescence vertical illuminator



Use of the fluorescent dyes

| Fluorochrome | Objects to be stained | Concentration and Solution | Staining time |
|-------------------------------|---|---|---|
| Diamino naphtylsulphonic acid | | Technique see R. C. NAIRN: Fluorescent Protein Tracing. E. & S. Livingstone Ltd., Edinburgh & London 3rd Edition 1969. Also contains further references. | |
| Berberin sulphate | Cell nuclei Bacteria | 1:1,000 1:1,000 | Aqueous Aqueous 10 sec 15–35 sec |
| Catecholamine | | For technique see FALCK B., N. A. HILLARP, G. THIEME and A. TORP: Fluorescence of catecholamines and related compounds condensed with formaldehyde. Journ. Histochem. and Cytochem. 10 (1962), 348. | |
| Euchrysin | Cell Nuclei Leucocytes Lymphocytes Mucus Wood | 1:10,000 1:1,000 1:20,000 1:1,000 | Aqueous Aqueous Aqueous Aqueous 3 min 4 min 5 min 15 sec |
| Primulin O | Wood Elementary virus bodies | 1:500 1:1,000 | Aqueous Aqueous 1 min 20 sec |
| Pyronin | | Techniques see H. C. BURCK: Histolog. Technik (1966), 123–124 Georg Thieme Verlag, Stuttgart. | |
| Quinacrine mustard | | For techniques see CASPERSSON T., L. ZECH, E. J. MODEST, G. E. FOLEY, U. WAGH and E. SIMONSSON: Chemical Differentiation with Fluorescent Alkylating Agents in VICIA FABA Metaphase Chromosomes. Experimental Cell Research 58 (1969), 128–140. | |
| Thioflavine S | Leucocytes Lymphocytes | 1:1,000 | Aqueous 4 min |
| Auroporphphin G | Mucus | 1:20,000 | Aqueous 5 min |
| Acridine yellow | Tubercle bacilli | 1:500 | Aqueous 30 sec |
| Acridine orange | Cell nuclei Mucus Vacuoles Bacteria | 1:10,000 1:20,000 1:5,000 1:30,000 | Aqueous Aqueous Alkaline Physiological saline solution 3 min 5 min 5–10 min |

Use of the fluorescent dyes

| Fluorochrome | Objects to be stained | Concentration and Solution | Staining time |
|---|---|--|---|
| Auramine | Tubercle bacilli | 1:1,000 Aqueous | 10 min |
| Coriphosphin O | Cell nuclei Leucocytes Lymphocytes Mucus Cell walls Wood Diphtheria bacteria | 1:10,000 1:1,000 1:20,000 1:1,000 1:1,000 1:1,000 Aqueous Aqueous Aqueous Aqueous | 3 min 4 min 5 min 1 min 1 min 10 min |
| Fluorescein isothiocyanate (FITC) | For technique see R. C. NAIRN: Fluorescent Protein Tracing. | | |
| Phosphin 3 R | Cell nuclei | 1:10,000 Aqueous | 3 min |
| Tetracycline | For technique see EGER W., F. GÖTZ and H. KÄMMERER: Herstellung von Dünnschliffen aus Knochen und Weichgeweben nach Markierung mit Tetracyclinen". Langenbecks Arch. Klin. Chir. 306 (1964), 205–214. | | |
| Lissamine-Rhodamine B (RB 200) | For technique see R. C. NAIRN: Fluorescent Protein Tracing. | | |
| Magdal red | Cell walls | 1:1,000 Aqueous | 1 min |
| Pararosaniline (Feulgen) | For technique see B. ROMEIS: Mikroskopische Technik 16th Edition (1968) 293–297. | | |
| Rhodamine B | Vacuoles | 1:1,000 Aqueous | |
| Acid Fuchsin | For technique see H. M. FROST: Bone remodelling dynamics. Ch. C. Thomas Springfield, USA-Illinois (1963). | | |
| Thiazin red R | Albumin (together with euchrysin) | 1:10,000 Aqueous | ½ min |
| Tetra-methyl rhodamine isothiocyanate (TRITC) | For technique see R. C. NAIRN: Fluorescent Protein Tracing. | | |

Fluorescence-free embedding media

| Embedding medium | Refractive index n_D |
|-------------------------------|------------------------|
| Entellan | 1.501 |
| Eukitt | 1.491 |
| Non-fluorescent immersion oil | 1.515 |
| Liquid paraffin | 1.480 |

Excitation and emission of the fluorescent dyes

| Fluorochrome | Exciting radiation | Mean exciting wave length (nm) | Mean wave length of fluorescence (nm) |
|--|--------------------|--------------------------------|---------------------------------------|
| Diamino naphthylsulphonic acid (DANS)* 5-hydroxy tryptamine (5-HT)* | Ultra-violet | 340 380–415 | 525 520–530 |
| Aurophosphin G | Violet/ blue | 450 | 580 |
| Berberin sulphate | | 430 | 550 |
| Catecholamine | | 410 | 470 |
| Euchrysin | | 430 | 540 |
| Primulin O | | 410 | 550 |
| Pyronin | | 410 | 540 |
| Quinacrine mustard* | | 430–460 | 490–530 |
| Thioflavine S | | 430 | 550 |
| Acridine yellow | blue | 470 | 550 |
| Acridine orange | | 470 | 530–650 |
| Auramine | | 460 | 550 |
| Coriphosphin O | | 460 | 575 |
| Fluorescein isothiocyanate (FITC)* | | 490 | 520 |
| Morin | | | |
| Phosphin 3 R | | 465 | 565 |
| Tetracycline | | 390 | 440 |
| Lissamine-Rhodamine B (RB 200)* | Green | 575 | 595 |
| Magdala red | | 540 | 570 |
| Pararosaniline (Feulgen) | | 570 | 625 |
| Rhodamine B | | 540 | 625 |
| Acid fuchsin | | 540 | 630 |
| Thiazin red R | | 510 | 580 |
| Tetramethyl rhodamine isothiocyanate (TRITC)* | | 540 | 570 |

To demonstrate primary fluorescence the specimen is first excited with ultra-violet light and observed through various suppression filters. If this does not produce results, continue examination with other exciting radiations according to the table.

The filter combination to be used in connection with the microscope outfits listed on pp. 28 often depends on the problem to be solved or the aim of the investigation. For this reason no fixed filter combinations for the individual fluorescent dyes are given here, but only the required exciting radiation.

For excitation one of the filter combinations associated with the required exciting radiation listed on p. 20 is used.

* Manufacturers' value. The other values were measured in stained specimens with the LEITZ microspectrograph.

Filter combinations (as per February 1972) +

| Exciting radiation | Broad-band transmitted light | Broad-band incident light | Narrow-band incident light | Selective excitation | Setting wheel on the illuminator | Dichroic beam splitting mirrors/built-in suppression filters | Suppression filters |
|--------------------|---|--------------------------------|----------------------------------|---|----------------------------------|--|---|
| Ultra-violet | 1mm UG 1 2mm UG 1 | 2mm UG 1 | 2x2mm UG 1 | | 1 | TK 400 / K 400 | from K 430 onwards corresponding to the fluorescence maximum |
| Violet | 3mm BG 3 | (2x) 3mm BG 3 | 3mm BG 3 + S 405 (Type AL) | KP 425 + 3mm BG 3 | 2 | TK 455 / K 460 | from K 460 onwards according to the fluorescence maximum |
| Blue | 1.5mm BG 12 3mm BG 12 2x3mm BG 12 | 3mm BG 12 | 1.5mm BG 12 + S 470 (Type AL) | FITC transmitted light K 490 | 3 | TK 510 / K 515 | from K 510 onwards according to the fluorescence maximum for FITC*: K 510 or S 525 Type AL |
| Green | | 2mm BG 36 + S 546 (Type AL) | KP 560 + 2mm BG 36 | 4 | TK 580 / K 580 | from K 580 onwards corresponding to the fluorescence maximum K 610 for TRITC* | |

With all filter combinations a 2mm KG 1 heat filter and a 4mm BG 38 red suppression filter must be inserted in the lamp housing. Filters inserted in the lamp housing are exposed to strong heat. It is therefore recommended to swing them out when they are not in use.

* With the two-wave-lengths excitation of FITC/TRITC double staining a suppression filter slide with the K 510 and K 610 or S 525 and K 610 filters side by side must be used.

** 2 KP 490 = KP 500

*** K 480 can be inserted in the lamp housing to eliminate autofluorescence excited by the UV portion of KP 490.



+ If the filter combinations or designations are changed an insert will be published for the filter table.

C Illuminating methods in fluorescence microscopy

Incident-light fluorescence

With fluorescence excitation with incident light the exciting radiation is directed at the specimen from above through an incident-light illuminator with dichroic beam splitting mirrors through the microscope objective. Here the full aperture of the objective is utilized for excitation. This method guarantees a particularly high fluorescence intensity because no light losses can occur through scattering or primary absorption in the specimen, which is possible with transmitted-light excitation.


Since the intensity of the exciting radiation depends on the numerical aperture of the objective used, microscope objectives of the largest possible aperture should be used for incident-light excitation.



Fluorescence vertical illuminator according to Ploem

- | | | | |
|---|---|----|----------------------------|
| 1 | Suppression filter slide | 5 | Light source |
| 2 | Suppression filter in the fluorescence illuminator | 6 | Filter in the lamp housing |
| 3 | Revolving turret with dichroic beam splitting mirrors and suppression filters | 7 | Field diaphragm |
| 4 | Dichroic beam splitting mirrors | 8 | Exciting filter turret |
| | | 9 | Stray-light stop |
| | | 10 | Objective |
| | | 11 | Specimen |

Later fluorescence vertical illuminators are equipped with an exciting filter turret (cf. illustration); this permits instant change of exciting filters.



Fluorescence vertical illuminator with exciting filter turret and filter holder with field diaphragm

For fluorometry incident-light illuminators with lateral mirror are used.

Built-in filters:

Position 1: 2mm UG 1

Position 2: 3mm BG 3 + S 405 interference band filter

Position 3: 2x KP 490 (= KP 500)

Position 4: 2mm BG 36 + S 546 interference band filter

Position 5: Empty aperture; use of any desired filter combinations possible.

Application

Investigations in biology and medicine (research and routine), investigation of opaque materials (natural sciences and technology).

It is unnecessary to centre the condenser, since the objective also fulfils the function of a condenser.

The objective aperture has the simultaneous function of field and aperture diaphragm.

Focusing on the specimen automatically produces quality of illumination (provided the lamp has been correctly centred).

Transmitted-light fluorescence


Fluorescence excitation with transmitted light is the classical method. It is carried out with darkground condensers or brightfield condensers.

Darkground condensers have the advantage that no exciting light reaches the objective and thus a dark background for the specimen is guaranteed. However, only objectives of apertures up to 1.15 can be used.

Brightfield condensers permit the use of objectives of larger aperture. But the image contrast is inferior, since here the exciting light also reaches the objective.

Simultaneous method

Simultaneous methods are used in fluorescence microscopy whenever besides the fluorescent structures the non-fluorescent parts of the object, too, are to be shown in brightfield, darkground, phase contrast or polarized light. The use of simultaneous methods is particularly favourable in connection with the fluorescence vertical illuminator, since here bright and contrasty fluorescence is preserved in the brightened field of view of the microscope.



27640-512 R

Fluorescence microscopy with darkground condensers

Application

Investigation in biology and medicine (routine and research)
Routine FITC fluorescence with KP 490 exciting filter and K 510 suppression filter.

Simultaneous incident-light and transmitted-light fluorescence.

Advantages over the brightfield condenser:

Dark background. Use of objectives with weak primary fluorescence possible.

Transmitted-light darkground / incident-light fluorescence


Applications

In biology and medicine

For the demonstration of the structure of the tissue surrounding the fluorescent parts of the object and for the localization of the fluorescent particles in tissue.

For the simultaneous demonstration of fluorescent and non-fluorescent object details whose colours can be also differentiated.

For the assessment of the number of non-fluorescent compared with that of the fluorescent particles.



27641-512 R

Transmitted-light darkground / transmitted-light darkground fluorescence


Applications

In biology and medicine

For the demonstration of the structure of tissue surrounding the fluorescent object parts and localization of the fluorescent particles in tissue.

For the simultaneous demonstration of fluorescent and non-fluorescent object details whose colours can also be differentiated.

For the assessment of the number of non-fluorescent compared with that of fluorescent particles.



Incident-light fluorescence / transmitted-light phase contrast


Applications

In biology and medicine

For the demonstration of the structure of the tissue surrounding the fluorescent object parts.

For the localization of the fluorescent particles in tissue.


For the assessment of the number of non-fluorescent compared with that of the fluorescent particles.



Incident-light fluorescence / transmitted polarized light

Applications


Simultaneous demonstrations of polarization and fluorescence in bi-refringent materials.



Simultaneous incident- and transmitted light fluorescence

Applications

Demonstration of two object portions with different fluorochromes (e.g. transmitted-light blue excitation and incident-light green excitation).



Incident-light fluorescence microscopes

| Microscope | Light source | Base Mirror housing | Filter* | Illuminator condenser | Objectives | Tube | Eyepieces |
|---|--|---|----------|---|--|------------------|-----------------------------------|
| DIALUX Routine | Lamp housing 100 Z with 50 W mercury lamp | Lamp holder | BB SA | Fluorescence vertical illuminator | Achromats, fluorite systems | S FSA | PERIPLAN® 6.3x PERIPLAN GF 10x |
| ORTHOLUX® ORTHOLUX II ORTHOPLAN Routine | Lamp housing 100 Z with 50 or 100 W mercury lamp; | ORTHOLUX: base A with mirror housing 30 S | BB SA | Fluorescence vertical illuminator according to PLOEM | Achromats, fluorite systems | FSA or FSA-GW | PERIPLAN 6.3x PERIPLAN GF 10x |
| | Lamp housing 250 with 200 W mercury lamp; | ORTHOLUX II: mirror housing 250/0 | | | | | |
| | ORTHOPLAN: lamp housing 500 with 200 W mercury lamp | ORTHOPLAN: mirror housing 250 S mirror housing 500 | | | | | |
| DIALUX ORTHOLUX ORTHOLUX II ORTHOPLAN Research | Lamp housing 100 Z with 100 W mercury lamp or 75 W xenon lamp | ORTHOLUX: base A with mirror housing 30 S | SB SA | Fluorescence vertical illuminator | Achromats, fluorite systems, apo chromats | FSA or FSA-GW | PERIPLAN 6.3x PERIPLAN GF 10x |
| | Lamp housing 250 with 200 W mercury lamp or 150 W xenon lamp | ORTHOLUX II: mirror housing 250/0 | | | | | |
| | ORTHOPLAN: Lamp housing 500 with 200 W mercury lamp 150 or 450 W xenon lamp | DIALUX: lamp holder | | | | | |
| ORTHOLUX ORTHOLUX II ORTHOPLAN Flurometry (MPV) | Lamp housing 100 Z with 100 W mercury lamp or 75 W xenon lamp | ORTHOLUX/ ORTHOLUX II: ARISTOPHOT®: ORTHOPLAN®: camera rail | SB SA | Fluorescence vertical illuminator according to PLOEM with lateral mirror | Achromats, fluorite systems, apo chromats | FS or FSA-GW | PERIPLAN 6.3x PERIPLAN GF 10x |

* BB = Broad-band excitation SB = Narrow-band excitation

SA = Selective excitation

Microscopes for fluorescence microscopy with darkground condenser

(Routine, research, fluorometry)

| Microscope | Light source | Base Mirror housing | Filter | Condenser | Objectives | Tube | Eyepieces |
|---|---|---|---|--|---|------------------|---|
| DIALUX® routine (FITC) | 12 v 100 W tungsten halogen lamp in lamp housing 100 or 100 Z | Feet with adapter | KP 490 + 4mm BG 38 / K 510 | Darkground condenser D 0.80 or D 1.20 | Achromats, fluorite systems | S FSA | PERIPLAN 6.3x PERIPLAN GF 10x |
| Routine (Broad band) | Lamp housing 100 Z with 50 W mercury lamp or lamp housing 250 with 200 W mercury lamp | Microscope base L or feet with adapter (only for lamp housing 100 Z) | BB | Darkground condenser D 0.80 or D 1.20 | Achromats, fluorite systems | S FSA | PERIPLAN 6.3x PERIPLAN GF 10x |
| ORTHOLUX ORTHOLUX II ORTHOPLAN Routine (FITC) | Lamp housing 100 Z with 50 W mercury lamp; lamp housing 250 with 200 W mercury lamp or 150 W xenon lamp ORTHOPLAN: lamp housing 500 with 150 or 450 W xenon lamp | ORTHOLUX: base A with mirror housing 30 S ORTHOPLAN: mirror housing 250 S | KP 490 + 4mm BG 38 / K 510 or S 525 | Darkground condenser D 0.80 or D 1.20 | Achromats, fluorite systems | FSA or FSA-GW | PERIPLAN 6.3x PERIPLAN GF 10x ORTHOPLAN PERIPLAN 6.3x GW/GG 8x, GW 10x |
| ORTHOLUX ORTHOLUX II ORTHOPLAN Routine (Broad band) | Lamp housing 100 Z with 50 W mercury lamp Lamp housing 250 with 200 W mercury lamp ORTHOPLAN: Lamp housing 500 with 200 W mercury lamp | ORTHOLUX: Base A with mirror housing 30 S ORTHOPLAN: mirror housing 250 S ORTHOPLAN: mirror housing 500 ORTHOLEX II: Mirror housing 250/0 | BB | Darkground condenser D 0.80 or D 1.20 | Apochromats, achromats fluorite systems | FSA or FSA-GW | PERIPLAN 6.3x PERIPLAN GF 10x ORTHOPLAN: PERIPLAN GW 6.3x GW/GG 8x, GW 10x |

Microscope for fluorescence microscopy with darkground condenser

(Routine, research, fluorometry)

| Microscope | Light source | Base Mirror housing | Filter | Condenser | Objectives | Tube | Eyepieces |
|----------------------|---|---|----------------|--|---|------------------|--|
| Research | Lamp housing 250 with 200 W mercury lamp or 150 W xenon lamp (FITC) ORTHOPLAN: Lamp housing 500 with 200 W mercury lamp, 150 or 450 W xenon lamp | ORTHOLUX: Base A with mirror housing 30 S ORTHOLUX II: Mirror housing 250/0 ORTHOPLAN: Mirror housing 250 S mirror housing 500 | BB SB SA | Darkground condenser D 0.80 or D 1.20 | Apochromats, fluorite systems, achromats | FSA or FSA-GW | PERIPLAN 6.3x PERIPLAN GF 10x ORTHOPLAN: PERIPLAN GW 6.3x GW/GG 8x, GW 10x |
| Fluorometry (MPV) | Lamp housing 250 with 150 W xenon lamp (if necessary for FITC lamp housing 100 Z with 12 v 100 W tungsten halogen lamp) | ORTHOLUX/ ORTHOLUX II: ARISTOPHOT: ORTHOPLAN: camera rail | SB SA | Darkground condenser D 0.80 or D 1.20 | Apochromats, fluorite systems | FS or FSA-GW | PERIPLAN 6.3x PERIPLAN GF 10x ORTHOPLAN: PERIPLAN GW 6.3x GW/GG 8x, GW 10x |

Fluorescence Microscopes for transmitted-light darkground/incident-light fluorescence

(Routine and research)

| Microscope | Light source | Base Mirror housing | Filter | Illuminator condenser | Objectives | Tube | Eyepieces |
|---|--|---|----------------|--|--|------------------|----------------------------------|
| DIALUX Routine | Lamp housing 100 Z with 50 W mercury lamp | Lamp holder | BB SA | Fluorescence vertical illuminator | Achromats, fluorite systems | S FSA | PERIPLAN 6.3x PERIPLAN GF 10x |
| ORTHOLUX ORTHOLUX II ORTHOPLAN Routine | Lamp housing 100 Z with 50 W mercury lamp Lamp housing 250 with 200 W mercury lamp or 150 W xenon lamp (FITC) ORTHOPLAN: Lamp housing 500 with 200 W mercury lamp | ORTHOLUX: base A with mirror housing 30 S ORTHOLUX II: Mirror housing 250/0 ORTHOPLAN Mirror housing 60 S (incident light) mirror housing 500 | BB | Fluorescence vertical illuminator according to PLOEM and darkground condenser D 0.80 or D 1.20 | Achromats, fluorite systems | FSA or FSA-GW | PERIPLAN 6.3x PERIPLAN GF 10x |
| DIALUX ORTHOLUX ORTHOLUX II ORTHOPLAN | Lamp housing 100 Z with 100 W mercury lamp or 75 W xenon lamp | ORTHOLUX: Base A with mirror housing 30 S | BB SB SA | Fluorescence vertical illuminator and dark- ground condenser | Achromats, fluorite systems Apochromats | FSA or FSA-GW | PERIPLAN 6.3x PERIPLAN GF 10x |
| Research | Lamp housing 250 with 200 W mercury lamp or 150 W xenon lamp ORTHOPLAN: Lamp housing 500 with 200 W mercury lamp, 150 or 450 W xenon lamp | ORTHOLUX II: Mirror housing 250/0 DIALUX: Lamp holder ORTHOPLAN: Mirror housing 60 S (incident light) Mirror housing 500 | | D 0.80 or D 1.20 | | | |

Fluorescence Microscopes for transmitted-light darkground/ transmitted-light darkground fluorescence

| Microscope | Light source | Base Mirror housing | Filter | Condenser | Objectives | Tube | Eyepieces |
|--|--|--|----------|---------------------------------------|--|---------------|--|
| DIALUX Routine | Lamp housing 100 Z with 50 W mercury lamp Lamp housing 250 with 200 W mercury lamp | Base L | BB | Darkground condenser D 0.80 or D 1.20 | Achromats, fluorite systems, | S | PERIPLAN 6.3x PERIPLAN GF 10x |
| ORTHOLUX II ORTHOLUX ORTHOPLAN Routine | Lamp housing 100 Z with 50 W mercury lamp Lamp housing 250 with 200 W mercury lamp ORTHOPLAN: Lamp housing 500 with 200 W mercury lamp | ORTHOLUX: Base A with mirror housing 30 S ORTHOLUX II: Mirror housing 250/0 ORTHOPLAN: Mirror housing 250 S ORTHOPLAN: Mirror housing 500 | BB | Darkground condenser D 0.80 or D 1.20 | Achromats, fluorite systems | FSA or FSA-GW | PERIPLAN 6.3x PERIPLAN GF 10x ORTHOPLAN: PERIPLAN GW 6.3x GW/GG 8x, GW 10x |
| Research | Lamp housing 250 with 200 W mercury lamp or 150 W xenon lamp (FITC) ORTHOPLAN: Lamp housing 500 with 200 W mercury lamp, 150 or 450 W xenon lamp | ORTHOLUX: Base A with mirror housing 30 S ORTHOLUX II: Mirror housing 250/0 ORTHOPLAN: Mirror housing 250 S ORTHOPLAN: Mirror housing 500 | BB SA | Darkground condenser D 0.80 or D 1.20 | Achromats, fluorite systems, apochromats | FSA or FSA-GW | PERIPLAN 6.3x PERIPLAN GF 10x ORTHOPLAN: PERIPLAN GW 6.3x GW/GG 8x, GW 10x |

Fluorescence Microscopes for incident-light fluorescence/transmitted-light phase contrast

| Microscope | Light source | Base Mirror housing | Filter | Condenser | Objective | Tube | Eyepieces |
|-------------|---|---|--------|--|-------------|---------------|---------------|
| DIALUX | Lamp housing | DIALUX: | BB | Fluorescence vertical illuminator and phase contrast condenser according to Heine or Zernike | Pv or PHACO | FSA or FSA-GW | PERIPLAN 6.3x |
| ORTHOLUX | 100 Z with | Lamp holder | | | | | |
| ORTHOLUX II | 50 W mercury | ORTHOLUX: | | | | | |
| ORTHOPLAN | lamp | Base A with mirror housing | | | | | |
| Routine | Lamp housing 250 with 200 W mercury lamp | 30 S | | | | | |
| | | ORTHOLUX II: Mirror housing 250/0 | | | | | |
| | ORTHOPLAN: Lamp housing 500 with 200 W mercury lamp | ORTHOPLAN: Mirror housing 60 S (incident light) Mirror housing 500 | | | | | |
| DIALUX | Lamp housing | DIALUX: | BB | Fluorescence vertical illuminator and phase contrast condenser according to Heine or Zernike | Pv or PHACO | FSA or FSA-GW | PERIPLAN 6.3x |
| ORTHOLUX | 100 Z with | Lamp holder | SB | | | | |
| ORTHOLUX II | 100 W mercury | ORTHOLUX: | SA | | | | |
| ORTHOPLAN | lamp or 75 W xenon lamp | Base A with mirror housing | | | | | |
| Routine | Lamp housing 250 with 200 W mercury lamp | 30 S | | | | | |
| | 150 W xenon lamp (FITC) | ORTHOPLAN: Mirror housing 250/0 | | | | | |
| | ORTHOPLAN: Lamp housing 500 with 200 W mercury lamp or 150 or 450 W xenon lamp | ORTHOLUX II: Mirror housing 60 S (incident light) Mirror housing 500 | | | | | |

Fluorescence Microscopes for incident-light fluorescence/transmitted polarized light

(Routine and Research)

| Microscope | Light source | Base Mirror housing | Filter | Condenser | Objectives | Tube | Eyepieces |
|---------------------|---|---------------------------------------|---|---|-------------------------------|------|----------------------------------|
| ORTHOLUX Routine | Lamp housing 100 Z with 50 W mercury lamp | Base A with mirror housing 30 S | BB and polarizer in holder and analyser in push-on mount | Fluorescence vertical illuminator and brightfield condenser series 600 | Polarized-light objectives | FSA | PERIPLAN 6.3x PERIPLAN GF 10x |
| Research | Lamp housing 100 Z with 100 W mercury lamp Lamp housing 250 with 200 W mercury lamp | Base A with mirror housing 30 S | BB SB and polarizer in holder and analyser in push-on mount | Fluorescence vertical illuminator and brightfield condenser series 600 (if necessary also polarized-light condenser) | Polarized-light objectives | FSA | PERIPLAN 6.3x PERIPLAN GF 10x |

Fluorescence Microscopes for simultaneous incident- and transmitted-light excitation

(Routine and Research)

| Microscope | Light source | Base Mirror Housing | Filter | Condenser | Objectives | Tube | Eyepieces |
|---------------------|---|---|----------------|--|--|------------------|----------------------------------|
| ORTHOLUX ROUTINE | Transmitted and incident light: Lamp housing 100 Z with 50 W mercury lamp Lamp housing 250 with 200 W mercury lamp ORTHOPLAN: for transmitted light: Lamp housing 250 on the mirror housing 500 with 200 W mercury lamp for incident light: 200 W mercury lamp in the lamp housing 500 | ORTHOLUX: Base A with 2 mirror housings 30 S For transmitted light: Mirror housing 60 ORTHOPLAN: in transmitted light: Mirror housing 250 S in incident light: Mirror housing 60 S | BB SA | Fluorescence vertical illuminator and darkground condenser D 0.80 or D 1.20 | Achromats, fluorite systems | FSA or FSA-GW | PERIPLAN 6.3x PERIPLAN GF 10x |
| Research | Incident light: Lamp housing 100 Z with 100 W mercury lamp or 75 W xenon lamp Transmitted light: Lamp housing 250 with 200 W mercury lamp or 150 W xenon lamp ORTHOPLAN: for transmitted light: Lamp housing 250 with 200 W mer- cury lamp or 150 W xenon lamp on mirror housing 500 for incident light: Lamp housing 500 with 200 W mercury lamp, 150 or 450 W xenon lamp | As under "Routine" | BB SB SA | Fluorescence vertical illuminator and darkground condenser D 0.80 or D 1.20 | Achromats, fluorite systems, apoachromats | FSA or FSA-GW | PERIPLAN 6.3x PERIPLAN GF 10x |

D Objectives

All objectives which because of their large aperture concentrate much exciting light on the specimen or produce a high yield of fluorescence light are particularly suitable. The intensities increase as the square of the apertures.

Example: With objectives of a given primary magnification but different apertures of 0.50 and 1.00 the fluorescent light yield of the objective of larger aperture is about 4 times as much as that of the smaller aperture objective.

Oil- or water immersion objectives are preferable. Here no loss of light owing to refraction on glass/air interfaces and no reflections between coverglass and front lens of the objective occur.*

For the simultaneous methods suitable phase contrast objectives or polarized light objectives must be used. Generally all objectives can be used in fluorescence microscopy. It must, however, be borne in mind that objectives of low magnification do not produce a very bright fluorescent image. They should therefore be used only for occasional low-power observation.

Transmitted-light fluorescence

| | |
|-------------------------|-------------------------|
| 10/0.25 | PI Apo 16/0.40 |
| 20/0.50 | Apo 25/0.65 |
| 40/0.65 | With adapter |
| 63/0.85 | PLEZY-FLU |
| With adapter | PI Apo 40/0.75 |
| PLEZY-FLU | Apo 63/0.95 |
| Oel 100/1.30 | with adapter |
| with funnel stop 513362 | PLEZY-FLU |
| | Oel 100/1.30 |
| | with funnel stop 513362 |

Incident-light fluorescence

| | |
|----------------------------|---------------------------------|
| Apo 25/0.65 | Oel + W 22/0.65 with |
| Fl 40/0.85 | adapter PLEZY-FLU*) |
| Apo 63/0.95 | Fl Oel 54/0.95 with |
| Fl Oel 95/1.32 | adapter PLEZY-FLU*) |
| or Iris Oel 100/1.30–1.10 | Pv Fl Oel 70/1.15 |
| | with adapter |
| | PLEZY-FLU*) |
| | Iris Fl Oel 95/1.32–1.10 |
| | or for QM chromosome |
| | fluorescence according |
| | to Caspersson |
| | Oel 100/1.30 |
| | with funnel stop 513362 |
| *) the adapter is required | only if objectives are used |
| | together with Oel 100 objective |
| | on the nosepiece. |

In addition the following water immersion objectives will shortly be available:

| |
|------------|
| W 25/0.60 |
| W 50/1.00 |
| W 100/1.20 |

They are particularly recommended for incident-light fluorescence.

The following objective combinations are particularly recommended:

E Eyepieces

The intensity of the fluorescent light depends on, among other factors, the magnification of the eyepieces. It decreases as the square of the eyepiece magnification increases.

Example: Eyepiece PERIPLAN 6.3x — intensity 100 %
Eyepiece PERIPLAN GF 12.5x — intensity 25 %.

Eyepieces of low magnification are therefore preferable for visual observation.

For fluorescence photomicrography (see Section III), however, the PERIPLAN GF 10x M or PERIPLAN GF 12.5x MF eyepieces with graticule for outlining the formats should be used for focusing.

III Fluorescence Photomicrography

A The specimen

should meet the following conditions:
Clean preparation without dust particles and air bubbles;
if possible no surplus dye in the interstices between fluorescent areas of the specimen (careful rinsing of the specimen during preparation, embedding in non-fluorescent medium);

Fluorescence as brilliant as possible:

Good differentiation of types of cell and detailed structures.

with full light intensity due to the sometimes low brightness.

For fluorescence photo micrography of weakly fluorescent specimens the monocular tube is particularly suitable: 100 % of the light is directed to the photo tube during photography, and 100 % to the eyepiece during observation.

Remove carefully dust from all glass surfaces in the beam path.

B The Microscope

The quality of the light source is important (before it is taken into use the lamp should be run for a day or longer; it will then burn more steadily with hardly any flicker).

When the intensity appreciably decreases or noticeable flickering occurs replace the lamp even before the end of its rated life.

Good centration of the light source and the condenser.

If possible use an objective of large aperture (immersion systems such as the FI Oel 54/0.95). The image brightness increases as the square of the objective aperture.

For visual observation choose eyepieces of the lowest possible magnification (image brightness decreases as the square of the eyepiece magnification increases). For fluorescence photomicrography eyepieces of high magnification and with graticules outlining the format are required (PERIPLAN GF 10x M, GF 12.5x MF).

If the FSA tube is used (beam 100 % or 20 % directed to the eyepiece, 80 % to the photographic device) it is recommended to focus the specimen in the eyepiece always

C The Camera

The following cameras are particularly well suited: The ORTHOMAT® automatic microscope camera with integrating and detail measurement for 35mm photography.

The large-format camera with the fully automatic exposure control for the 9 x 12cm (4 x 5in) format.

For brightly fluorescent specimens the micro-attachment for the LEICA®, the system attachment camera and the camera with POLAROID* CB 100 back are suitable; precise exposure measurement is carried out with the MICROSIX®-L exposure meter. If its measuring range is inadequate, an exposure series should be taken, if necessary.

For the accurate determination of the exposure time the measuring field of the exposure meter or of the camera should be completely occupied by a specific portion of the specimen.

If there is a tendency of underexposure because of the bleaching-out of small particles set the object ratio of the measuring field (e.g. ORTHOMAT) at 80–100 %.

* = Trade mark of Polaroid GmbH.

D The Exposure Material

For fluorescence photomicrography daylight colour films are to be preferred to black-and-white films because of their reproduction of the specific colour of fluorescence*.

The results can be compared and reproduced because of the standardized development of colour reversal films.

Reproduction from colour originals is also recommended for black-and-white reproduction in print. Colour enlargements can also be obtained from colour transparencies.

Colour negative films are not recommended because of their generally low sensitivity and the subjectivity of positive production. Usually processing is carried out by processing stations where the technicians are unfamiliar with the fluorescence phenomenon. If these films are used nevertheless, a suitable transparency should, if possible, be included with the film when it is sent to the processing station to make colour comparisons possible.

Black-and-white films should be developed for good contrast.

The following exposure material was used in our Applied Microscopy Laboratory:

Black-and-white films

| Manufacturer | Type of film | Speed (DIN) (ASA) | Remarks |
|--------------|-----------------------|---------------------------|-----------------------------------|
| Agfa-Gevaert | Scientia 50 B 65 | 18 50 | |
| Kodak | Plus-X Tri-X | 22 125 27 400 | very good contrast and fine grain |
| Ilford | Pan F FP 4 HP 4 | 18 50 22 125 27 400 | |

Colour films

| | | | |
|--------------|------------------------------|--------|--|
| Agfa-Gevaert | Agfachrome professional 50 S | 18 50 | daylight film |
| | Agfachrome professional 50 L | 20 80 | artificial-light film; optimum blue rendering (e.g. catecholamine fluorescence) |
| Kodak | CT 18 | 18 50 | daylight film |
| | Ektachrome High Speed | 23 160 | daylight film; special development when the speed was set at 27 DIN (400 ASA) |
| GAF | Ansco-chrome 500 | 28 500 | daylight film; recommended for the photography of extremely small objects (e.g. spirochaetes), since the fluorescence intensity quickly decreases owing to bleaching |
| ORWO | UT 21 | 21 100 | daylight film |

* For the photography of blue fluorescent colours (e.g. catecholamine/neurotransmitter fluorescence), however, artificial-light colour films should be used because of their superior blue rendering.

IV References

- BECKER, H.: „Fluoreszenzhistochemie biogener Amine“. Mikroskopie Bd. 26 (1970), S. 153–168.
- BERTALANFFY, L. v. and BERTALANFFY, F. D.: „Evaluation of the acridine-orange fluorescence microscope method for cytodiagnosis of cancer.“ Annals of the New York Academy of Sciences, Vol. 93, Art. 16, 715–750, Sept. 1962.
- BÖHM, N. und E. SPRENGER: „Fluorescence Cytophotometry: A Valuable Method for the Quantitative Determination of Nuclear Feulgen DNA.“ Histochemistry 16 (2), 100–118, 1968.
- BRÄUTIGAM, Fr. und GRABNER, A.: „Beiträge zur Fluoreszenz-Mikroskopie.“ Georg Fromme, Wien, 1949.
- CASPERSSON, T., ZECH, L., MODEST, E. J., FOLEY, G. E., WAGH, U. und SIMONSSON, E.: „Chemical Differentiation With Fluorescent Alkylation Agents in VICA FABA Metaphase Chromosomes.“ Experimental Cell Research, 58 (1969), S. 128–140.
- CASPERSSON, T., ZECH, L., JOHANSSON, C. und MODEST, E. J.: „Identification of Human Chromosomes by DNA-Binding Fluorescent Agents.“ Chromosoma (Berl.) 30 (1970), S. 215–227.
- CHADWIG, C. S., McENTEGART, M. G. and NAIRN, R. C.: (a): „Fluorescent Protein Tracers, a simple alternative to fluorescein“; Lancet 1 (1958), 412–414.
(b): „Fluorescent Protein Tracers, a trial of new fluorochromes and the development of an alternative to fluorescein“; Immunology 1 (1958), 315–327.
- CHERRY, W. B., GOLDMAN, M. and CARSKI, T. R.: „Fluorescent Antibody Techniques in the diagnosis of Communicable Diseases“; U. S. Dept. of Health, Education, and Welfare, Publ. Health Serv. Publication No. 729, 1960.
- CLAYTON, R. M.: „Localization of embryonic antigens by antisera labelled with fluorescent dyes“; Nature (Lond.) 174 (1954), 1059.
- CONTIER, L.: „Die Vitalfluorochromierung der Plasmazytomzellen mit Acridinorange“; Blut 3, 185–191 (1957).
- COONS, A. H. and KAPLAN, M. H.: „Localization of antigen in tissue cells. II. Improvements in a method for the detection of antigen by means of fluorescent antibody.“ The Journal of Experimental Medicine, II, 1–13 (1950).
- COONS, A. H.: „Fluorescent antibodies as histochemical tools.“ Federation Proceedings 10, Nr. 2, 558–559 (1951).
- CORRODI, H. and J. JONSSON: „The Formaldehyde Fluorescence Method for the Histochemical Demonstration of Biogene Monoamines.“ Journ. Histochem. and Cytochem. 15 (1967), 65.
- DANIELLI, J. F.: „General cytochemical methods.“ Academic Press, New York. 1958.
- DE LERMÄ, B.: „Die Anwendung von Fluoreszenzlicht in der Histochemie.“ S. 78–159 (ca. 1000 Literaturzitate im Anhang!) in Handbuch der Histochemie I/1. G. Fischer Verlag. Stuttgart. 1958.
- EDER, H.: „Fluoreszenzmikroskopische Eosinophilen-Diagnostik.“ LEITZ-Mitt. Wiss. u. Techn. I, (1960) 113.
- EDER, H. und FRITSCHE, H.: „Zur Farbphotographie fluoreszenzmikroskopischer Objekte“; LEITZ-Mitt. Wiss. u. Techn. II (1963), 143–145.
- EICHLER, J. und F. WALTER: „Ein Beitrag zur Fluoreszenzmikroskopie des Knochengewebes“; LEITZ-Mitt. Wiss. u. Techn. IV (1967), 110.
- ERÄNKÖ, O.: „Distribution of fluorescing islets, adrenaline and noradrenaline in the adrenal medulla of the hamster.“ Acta endocrinol. 18 (1955), 174.
- FALCK, B., N. A. HILLARP, G. THIEME and A. TORP: „Fluorescence of catecholamines and related compounds condensed with formaldehyde“; Journ. Histochem. and Cytochem. 10 (1962), 348.
- GREHN, J. and H. KORNMANN: „Kontrastfluoreszenz mit Opak-Illuminator“; LEITZ-Mitt. Wiss. u. Techn. III (1965), 108.
- GOTTSCHEWSKI, G. H. M.: „Apparate und Einrichtungen für qualitative fluoreszenzmikroskopische Untersuchungen“; S. 160 bis 170 in Handbuch der Histochemie I/1. G. Fischer Verlag. Stuttgart. 1958.
- HAGEMANN, P. K. H.: „Fluoreszenzfärbung von Tuberkelbakterien mit Auramin“; Münchener Medizinische Wochenschrift 85, 1066–1068 (1938).
- HAITINGER: „Fluoreszenz-Mikroskopie.“ 2. Auflage. Akademie-Verlag. Leipzig. 1960.

- HAMPERL, H.: „Microscopia con luz fluorescente.“ Trab. Inst. Cajal invest. biol. 48, 253–271 (1956).
- HANSEN, P. A.: „Fluorescent Compounds used in Protein Tracing Absorption and Emission Data.“ Arbeitsbericht des Dept. of Microbiology, University of Maryland. Sept. 1964 (nicht im Handel erhältlich).
- HARDERS-STEINHÄUSER, M. und BAUER, G.: „Bewährung der vollautomatischen Mikroskopkamera ORTHOMAT bei der Fluoreszenzmikroskopie“; LEITZ-Mitt. Wiss. u. Techn. II (1963), 182–183.
- HIJMANS, W., H. R. E. SCHUIT, A. P. M. JONGSMA and J. S. PLOEM: „Performance Testing of Fluorescent Antisera Against Human Immunoglobulins“, 1970 in E. J. HOLBOROW (Ed.) Standardization in Immunofluorescence 193–202, Blackwell Scientific Publications, Oxford and Edinburgh 1970.
- KING, N.: „Fluorescence microscopy of fat in milk and milk powder.“ The Journal of Dairy Research 22, 205–210 (1955).
- KING, J.: „Fluorescence Microscopy.“ S. 34–42 in A. E. J. Vickers, Modern methods of microscopy. Butterworths Scientific Publications (1956).
- KOCH, E. und KONOFSKY, B.: „Die Eosinophilenzählung mit der Fluoreszenzmethode nach Eder.“ Die Medizinische Welt Nr. 14, 721–723 (1960).
- KOCH, K.-F.: „Lichtquellen für die Fluoreszenzmikroskopie. 1. FITC-Immunfluoreszenz.“ LEITZ-Mitt. Wiss. u. Techn. V (1971), 146–148.
- KOSENOW, W. und E. SCHELLOW: „Gleichzeitige Kammerzählung und -differenzierung der Leukozyten im Fluoreszenz-Mikroskop.“ Blut 4, 66–77 (1958).
- KOSENOW, R.: „Lebende Blutzellen im Fluoreszenz- und Phasenkontrastmikroskop“; S. Karger. Basel, New York. 1956.
- KRAFT, W.: „Die Technologie des LEITZ-Fluoreszenzopak“; LEITZ-Mitt. Wiss. u. Techn. IV (1969), 239.
- KRAFT, W.: Ein neues FITC-Erregerfilter für die Routinefluoreszenz. LEITZ-Mitt. Wiss. u. Techn. V, (1970), 41–44.
- KRIEG, A.: „Fluoreszenzanalyse und Fluorochromie in Biologie und Medizin“; Klinische Wochenschrift 31 (15), 350–356 (1953).
- MAIER, D. und WETZLE, W.: „Fluoreszenz-Mikroskopie geologischer und paläontologischer Objekte“; Zeiss-Mitt. I (4), 127–131 (1958).
- MAYERSBACH, H. v.: „Immunhistologische Methoden in der Histochemie“; S. 599–639 in Handbuch der Histochemie I/1. G. Fischer Verlag. Stuttgart. 1958.
Dort weitere Literatur.
- MEISSEL, M. N., KABANOVA, E. A., LEVINA, E. N. and STRACHOWA, V. A.: „Some new ways of fluorescence microscopy in the microbiology.“ Izv. Akad. Nau SSSR, Ser. Biol. 1958, 533–543 (russ. mit engl. Zusammenfassung).
- NAIRN, R. C.: „Fluorescent Protein Tracing.“ E. & S. Livingstone Ltd.; Edinburgh and London. 1. Aufl. 1962; 3. Aufl. 1969.
Dort weitere Literatur.
- NAIRN, R. C., CHADWICK, C. S. and McENTEGART, M. G.: „Fluorescent protein tracers in the study of experimental liver damage“; Journ. Path. Bact. 76 (1958), 143–153.
- NITSCH, B., MURKEN, J.-D. und BRÜCK, H.-J.: „Determinating Feulgen-DNA of Individual Chromosomes by Fluorescence Cytophotometry with Incident Light.“ Histochemistry 23 (1970), S. 254 bis 265.
- ÖYSTÈSE, B.: „Zur Kritik der fluoreszenz-immunologischen Methoden in der Histologie.“ Verh. dtsch. Ges. Pathol. (1962), 120 bis 125 (46. Tagung der Ges.).
- PEARSON, P. L. und BOBROW, M.: „Technique of Identifying Y-Chromosomes in Human Interphase Nuclei.“ Nature 226 (1970), April 4.
- PERNER, E.: „Die Methoden der Fluoreszenz-Mikroskopie“ in Freund, H., Handbuch der Mikroskopie in der Technik I/1, S. 357 bis 431.
- PERNIS, B., FORNI, L. und AMATE, L.: „Immunoglobulin spots on the surface of rabbit lymphocytes.“ J. exp. med. 132 (1970).
- PETERS, Th.: „Apparatur und Technik zur Mikroskopie an lebenden Säugeorganen in situ in gewöhnlichem Licht und Fluoreszenzlicht.“ Z. wiss. Mikrosk. 62, 348–368 (1955).
- PLOEM, J. S.: „A new microscopic method for the visualisation of blue formaldehydeinduced catecholamine Fluorescence“; Arch. int. Pharmacodyn. 182/2 (1969), 421.
- PLOEM, J. S.: „The Use of a Vertical Illuminator with Interchangeable Dichroic Mirrors for Fluorescence Microscope with incident Light.“ Z. Wiss. Mikroskopie 68 (3), 129–142 (1967).

PLOEM, J. S.: „Quantitative Immunofluorescence“, 63–73. Standardization in Immunofluorescence in E. J. HOLBOROW (Ed.), Blackwell Scientific Publications, Oxford and Edinburgh 1970.

PLOEM, J. S.: „Standards for Fluorescence Microscopy“, 137–153. Standardization in Immunofluorescence in E. J. HOLBOROW (Ed.), Blackwell Scientific Publications, Oxford and Edinburgh 1970.

PLOEM, J. S.: „Ein neuer Illuminator-Typ für die Auflicht-Fluoreszenzmikroskopie“; LEITZ-Mitt. Wiss. u. Techn. IV (1969), 225.

PLOEM, J. S.: „Die Möglichkeiten der Auflichtfluoreszenzmethode bei Untersuchungen von Zellen in Durchströmungskammern und Leightonröhren.“ 10. Symp. d. Gesellsch. f. Histochemie, Nijmegen 1965, Acta Histochem. Suppl. 7, 339–343 (1967).

PRICE, G. R. und S. SCHWARTZ: „Fluorescence microscopy“, in OSTER, G. and A. W. POLLISTER: Physical Techniques in biological Research. Vol. III, S. 91–139. Academic Press. N. Y. 1956.

RICHARDS, O. W. and P. WATERS: „A new Interference Exiter Filter for Fluorescence Microscopy of Fluorescein-Tagged Substances.“ Stain Technology 42 (1967), 320.

RIGGS, J. L. et al.: „Isothiocyanate compounds as fluorescent labelling agents for immune serum“; Amer. Journ. Path. 34 (1958), 1081–1097.

RYGAARD, J. and W. OLSEN: „Interference Filters for improved Immunofluorescence Microscopy.“ Acta path. microbiol. scand. 76 (1969), 146.

SCHRAMM, G. und RÖTTGER, B.: „Untersuchungen über das Tabakmosaikvirus mit fluoreszierenden Antikörpern“; Zeitschr. f. Naturforschung 14b (8/9), 510–515 (1959).

SCHUIT, H. R. E.: „Photomicrographic Recording of Immunofluorescence.“ In E. J. HOLBOROW (Ed.), Standardization in Immunofluorescence (1970), 159–162. Blackwell Scientific Publications, Oxford and Edinburgh 1970.

SCHÜMMELFEDER, N.: „Fluoreszenzmikroskopische Histochemie und Farbphotographie“; Photographie und Wissenschaft 6, 3/4, 19–24 (1957).

SERNETZ, M., and A. THAER: „A capillary fluorescence standard for microfluoreometry.“ Journal of Microscopy 91, 43–52, 1970.

STOCKINGER, L.: „Fluoreszenzuntersuchungen an Gewebekulturen. I. Die metachromatische Vitalfluoreszenz mit Acridinorange.“ Z. Naturforschg. 13b, 407–409 (1958).

STRUGGER, S.: „Fluoreszenz-Mikroskopie und Mikrobiologie“; Schaper, Hannover. 1949.

TCHAN, Y. T.: „Phase-fluorescence microscopy.“ Nature 179, 1243 bis 1244 (1957).

THAER, A.: „Microfluometric Analysis of the Reticulocyte Population in Peripheral Blood of Mammals.“ Cytol. Automation, Proceedings of the sec. Tenovus Symposium. Cardiff 1968, 189 bis 195. London/Edinburgh (Livingstone) 1970.

WALTER, F.: „Über die Fluoreszenzmikroskopie mit markierten Proteinen.“ LEITZ-Mitt. Wiss. u. Techn. II (1964), 207.

WALTER, F.: „Eine Auflicht-Fluoreszenz-Einrichtung für die Routinediagnose.“ LEITZ-Mitt. Wiss. u. Techn. IV (1968), 186.

WALTER, F.: Fluoreszenzmikroskopie in Biologie und Medizin. LEITZ-Mitt. Wiss. u. Techn. V, (1970), 33–40.

WALTER, F. und EICHLER, J.: „A new fluorescence microscope technique applied to bone tissues.“ Vortrag anl. Tagung der „Royal Microscopical Soc. and Anatomical Soc. of Great Britain and Ireland“ in Sheffield (England), April 1968.

WITTEKIND, D.: „Über die derzeitige Bedeutung der Vitalfärbung, speziell der Vitalfluorochromierung in der Cytologie und über die Möglichkeiten ihrer weiteren Anwendung, insbesondere in Kombination mit anderen Methoden.“ Mikroskopie 14, 9–25 (1959).

WUNDERLICH, H.: „Blutzählkammern mit fluoreszierender Netzteilung für Fluoreszenzuntersuchungen“; Z. wiss. Mikroskop. 64, 47–49 (1958).


YOUNG, M. R. and J. A. ARMSTRONG: „Fluorescence Microscopy with the Quartz-Jodine-Lamp.“ Nature 213 (1967), 649.

V List of LEITZ literature on the instruments mentioned in this booklet


| | | |
|--|--------------|---------|
| Image-forming and Illuminating Systems of the Microscope | List | 512-99 |
| Incident-light fluorescence | List | 512-105 |
| 4 x 5in bellows camera with fully automatic exposure control | List | 540-28 |
| | Instructions | 540-29 |
| DIALUX | List | 512-111 |
| | Instructions | 512-112 |
| The Microscope and its application | List | 512-69 |
| Fluorescence vertical illuminator according to Ploem | Instructions | 512-92 |
| 200 W Ultra-high pressure mercury lamp in the lamp housing 250 | Instructions | 514-117 |
| Lamp housing 250 | Instructions | 514-72 |
| Lamp housing 100/100 Z | Instructions | 514-119 |
| Micro attachment for the LEICA with vibration damper | List | 540-22 |
| | Instructions | 54-24 |
| MICROSIX-L exposure meter | List | 540-22 |
| | Instructions | 540-21 |
| ORTHOLUX | List | 512-40 |
| | Instructions | 512-76 |
| ORTHOMAT | List | 540-20 |
| | Instructions | 540-19 |
| ORTHOPLAN | List | 512-82 |
| | Instructions | 512-83 |
| Phase contrast devices | List | 513-5 |
| Phase contrast equipment with condenser according to Heine | Instructions | 513-21 |
| Phase contrast equipment with condenser according to Zernike | Instructions | 513-84 |
| Phase contrast fluorescence device | Instructions | 513-85 |

VI Index

| | |
|--|------------|
| Artificial-light colour film | 39 |
| Basic stands for fluorescence | 6 |
| Black-and-white films | 39 |
| Brightfield condenser | 23 |
| Broad-band excitation | 10 |
| Cameras | 38 |
| Colour negative film | 39 |
| Cytochemical staining | 5 |
| Darkground condenser | 23 |
| Daylight colour film | 39 |
| Dichroic beam splitting mirrors | 16, 20 |
| Excitation | 19 |
| Exciting filters | 13, 20 |
| Exposure material | 39 |
| Eyepieces | 37 |
| Filters | 10 |
| Filter combinations | 20 |
| Filter curves | 13 |
| Fluorescence | 4 |
| Fluorescence vertical illuminator according to Ploem | 21 |
| Fluorescent dyes | 17, 19 |
| Fluorochroming | 5 |
| Heat filter | 15, 20 |
| Histochemical staining | 5 |
| Immunofluorescence | 5 |
| Incident-light fluorescence | 21 |
| Instruments | 6, 29 |
| Intensity comparison of light sources | 8 |
| Interference band filters | 11, 13, 14 |
| Interference selection filters | 12, 16 |
| Light sources | 8 |
| Literature | 40 |
| Luminescence | 4 |
| Microscope | 7, 38 |
| Monocular tube | 38 |
| Narrow-band excitation | 11 |
| Objectives | 36 |
| Objective aperture | 36 |
| Oil immersion objectives | 36 |
| ORTHOMAT | 38 |
| Phosphorescence | 4 |
| Photomicrography | 38 |
| Primary filters | 7 |
| Primary fluorescence | 5 |
| Principle of construction of the fluorescence microscope | 7 |
| Red suppression filter | 12, 15 |
| References | 43 |
| Secondary filters | 7 |
| Secondary fluorescence | 5 |
| Selective excitation | 12 |
| Simultaneous method | 23 |
| Specimen | 38 |
| Spectral energy distribution HBO 200 | 8 |
| Suppression filters | 15 |
| Technical data of light sources | 9 |
| Transmitted-light fluorescence | 23 |
| Tungsten halogen lamp | 9 |
| Ultra-high pressure mercury lamps | 9 |
| Water immersion objectives | 36 |
| Xenon lamps | 9 |



1
Muscle, snail
Catecholamine
Fl Oel 54/0.95
216x
Incident-light
fluorescence
BG 3 + S 405, K 460




4
Eosinophilia
Aniline blue/
acridine orange
Oel 100/1.30
with funnel stop
400x
UG 1, K 430
Incident-light
fluorescence




2
Muscle, snail
Catecholamine
+ 5-HT, otherwise
as Fig. 1




5
Streptococci
Group B (negative
demonstration)
Apo Oel 90/1.40
360x
Incident-light
fluorescence
BG 12 + S 470/K 510




3
Mesentery, rat
Catecholamine
PI 16/0.40
64x
Incident-light
fluorescence
BG 3 + S 405, K 460




6
Tooth, Pomeranian
Tetracycline
PI Fl 4/0.14
16x
Transmitted-light
fluorescence
BG 12, K 530




7
Tooth, Pomeranian
Tetracycline,
otherwise as Fig. 6




10
Intestine,
guinea pig
Acridine orange
Oel + W 22/0.65
88x
Transmitted-light
fluorescence
BG 12, K 530




8
Kidney, rat
Acridine orange
25/0.50
100x
Transmitted-light
fluorescence
BG 12, K 530




11
Testicle, rat
Acridine orange
Fl Oel 54/0.95
216x
Incident-light
fluorescence
BG 12, K 530




9
Bacteria on
Epithelial cells of
the buccal mucosa
Acridine orange
40/0.65
160x
Transmitted-light
fluorescence
BG 12, K 530




12
FITC demonstration
of measles
Oel + W 22/0.65
88x
Incident-light
fluorescence
2x KP 490, S 525



13
Streptococci
Group B (positive
demonstration)
Apo Oel 90/1.40
360x
Incident-light
fluorescence
BG 12 + S 470/K 510




14
FITC demonstration
of German measles
Fl Oel 54/0.95
216x
Transmitted-light
fluorescence
KP 490, K 510




15
Gonococci (positive
demonstration)
in pus cell
Oel 100/1.30
400x
Incident-light
fluorescence
BG 12 + S 470/K 510


16
FTA-test
FITC
Apo Oel 90/1.40
360x
Incident-light
fluorescence
2x KP 490, K 510




17
FTA-test
FITC
Fl Oel 54/0.95
216x
Incident-light
fluorescence
2x KP 490, K 510




18
Bone marrow cells
FITC
Fl Oel 54/0.95
216x
Incident-light
fluorescence
2x KP 490, S 525







19
Bone marrow cells
FITC
Fl Oel 54/0.95
216x
Incident-light
fluorescence
2x KP 490, S 525




22
FITC demonstration
of measles
Oel + W 22/0.65
88x
Transmitted-light
fluorescence
BG 12, S 525




20
Kidney glomerulus
FITC
25/0.50
100x
Transmitted-light
fluorescence
BG 12, K 510



23
Bone marrow cells
(double labelling
with FITC and
TRITC)
Fl Oel 54/0.95
216x
Right:
2 KP 490/K 510
Centre:
2 KP 490/S 525
Left:
S 546 + BG 36/K 610



21
Bacteria coli
FITC
Fl Oel 95/1.32
380x
Incident-light
fluorescence
BG 12, K 510



24
Onion root
Feulgen
Fl Oel 54/0.95
216x
Transmitted-light
fluorescence
S 546 + BG 36, K 610



MICRO SCIENTIFIC CO., INC.

P. O. BOX 476
5701 HIGH ST.
WORTHINGTON, OHIO 43085
614 - 888 - 5697

® = Registered Trademark

Design subject to alterations without notice.

ERNST LEITZ GMBH D-6330 WETZLAR GERMANY

Subsidiary: Ernst Leitz (Canada) Ltd., Midland, Ontario