

The goniometer eyepiece (15× compensated eyepiece) is used for angular measurements. Instead of a micrometer, the glass plate behind the collector lens carries a circle, divided into 360°. The eyelens mount can be turned, whereas the one of the collector lens can be fixed at given rotation angles by means of a clamping screw on the tube. When turning the eyelens, a reticle cross in the field of view rotates over the fixed scale of angles. Reading is to one half of a degree, while fifths can still be estimated.


15× goniometer eyepiece (fig. 7)

Key number
5670



Fig. 7. Goniometer eyepiece.

Fig. 8. Field of view of the 15× goniometer eyepiece.



Illustrations and descriptions contained herein are not binding for sales contracts. Printed in Switzerland - M1 623 e - IX. 60

Accessories for microscopic measurements

Since the use of microscopy for technical purposes is becoming more and more frequent, precise microscopic measuring devices prove to be of ever increasing importance in almost any field of industrial enterprise. In the field of microscopic research also, precise measurements are often needed.

Our manufacturing program comprises the following special objectives, micrometer plates and object micrometers for measurements of lengths, angles and surfaces:

Object micrometer

Key number
8410

Scale of 5 mm divided into 1/2 mm, 2 mm divided into 1/10 mm, 0,2 mm divided into 1/100 mm, in case (fig. 1)
This object micrometer serves for calibrating purposes, e.g. for the determination of the length to which corresponds one part of the division of the micrometer plate in the object (as to the micrometer value, see instructions "microscopic measuring methods").



Fig. 1. Object micrometer.

Huygens type measuring eyepieces without micrometer plate

6× measuring eyepiece (field number 18) 5652
10× measuring eyepiece (field number 14) 5655

Huygens type measuring eyepieces with micrometer plate

6× measuring eyepiece with 10:100 plate (10 mm divided into 100 parts, fig. 2) 5653
10× measuring eyepiece with 5:100 plate (5 mm divided into 100 parts, analogous to fig. 2) 5656

Fig. 2. Field of view of the single eyepiece measuring plate division 5:100.





Fig. 5. Eyepiece with micrometer screw.

Eyepiece with micrometer screw

Key number
5675


For length measurements of highest accuracy the 15× eyepiece with micrometer screw is used (compensated eyepiece, 6 mm divided into 12 parts)

The approximate micrometer values of the measuring plates as listed above, resulting with various objectives, are given in the instruction booklet "microscopic measuring methods".

Single eyepiece measuring plates (Ø 16 mm)

Key number
5680
5681
5682
5683
5684

- Scale 10 mm divided into 100 parts (fig. 2)
- Scale 5 mm divided into 100 parts (analogous to fig. 2)
- Rectangular division (5 mm²) divided into squares of 1 mm side length
- Rectangular division (10 mm²) divided into squares of 1 mm side length (fig. 3)
- Cross hairs (fig. 4)



Field of view of the single eyepiece measuring plate.
Fig. 3. Rectangular division (10 mm²).
Fig. 4. Cross hairs.




Fig. 4




Fig. 6. Field of view of the eyepiece with micrometer screw. The cross hairs can be moved with respect to the field of view by actuating the measuring drum.